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URSL - 0400	2 x 62,5 mL
URSL - 0420	4 x 62,5 mL
URSL - 0500	5 x 125 mL
URSL - 0407	2 x 62,5 mL
URSL - 0427	4 x 62,5 mL
URSL - 0507	5 x 125 mL

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English - GB

In vitro diagnostic reagent, for professional use only

CLINICAL SIGNIFICANCE (1-3)

Urea is the major metabolite product of protein catabolism. The biosynthesis of urea from ammonia is exclusively carried out by hepatic enzymes. More than 90% of urea is excreted through the kidneys, with the remainder excreted through the gastrointestinal tract or skin. Blood urea concentrations can be increased by numerous factors linked to pre-renal causes (increased protein catabolism, as in haemorrhage into gastrointestinal tract, shock, some chronic liver diseases) or renal/postrenal causes (acute or chronic renal diseases, postrenal obstruction to urine flow). Uremia is also increased by high-protein diet, state of dehydration, muscle wasting (as in starvation). The determination of urea rate is used together with the determination of creatinine rate to discriminate between pre-renal (normal creatinin) and renal/postrenal (increased creatinin) disorders.

METHOD (3-5)

Enzymatic - UV.
Kinetic.

PRINCIPLE (4)

Enzymatic determination of urea according to the following reactions :

L-Glu = L-Glutamate
GIDH = Glutamate dehydrogenase

REAGENTS COMPOSITION

Reagent : R	
Tris buffer, pH 7,60 (37 °C)	125 mmol/L
ADP	1 mmol/L
α -Ketoglutarate	9 mmol/L
Uréase	≥ 8100 U/L
GIDH	≥ 1350 U/L
Reagent : R2	
NADH	1,5 mmol/L
Standard: Std	
Urea	50 mg/dL
	8,33 mmol/L

This standard is included in kits URSL-0407/0427/0507 and can be sold separately under the reference URUV-0055.

MATERIAL REQUIRED BUT NOT PROVIDED

CONT-0060	ELITROL I	10 x 5 mL
CONT-0160	ELITROL II	10 x 5 mL

PRECAUTIONS

- The reagents contain less than 0.1% of sodium azide. Sodium azide can react with copper and lead plumbing to form explosive metal azides. If discharge in the canalisations, rinse with plenty of water.

- Use clean or single use laboratory equipment only to avoid contamination.

- The standard should be immediately and tightly capped to prevent contamination and evaporation.

- For more information, Material Safety Data Sheet (MSDS) is available on request for professional user.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

STABILITY OF REAGENTS

Store at 2-8 °C and protect from light.

The reagents and the standard are stable until the expiry date stated on the label.

On board stability :

The stability is specific for each analyser (for COBAS MIRA refer to § PERFORMANCE DATA)

PREPARATION AND STABILITY OF WORKING REAGENT

The standard is ready to use.

- **One-reagent procedure**

Mix 4 volumes of reagent R1 with 1 volume of reagent R2.

Stability : 5 days at 20-25 °C
4 weeks at 2-8 °C

- **Two-reagent procedure**

The reagents are ready to use.

SAMPLES (2,3)

- **Specimen**

- Serum or heparinized plasma (except ammonium heparin). Do not use fluoride as inhibitor of glycolysis since it inhibits urease.

- Urines should be diluted at 1/20 to 1/50 with distilled water before analysis (when there is no pre-dilution by the analyser).

- **Storage**

Serum and plasma are stable up to 24 hours at room temperature, for one week at 4 °C. Frozen between -15 and -20 °C, these samples are stable for at least 2-3 months.

Urine samples are stable up to 4 days if stored at 4-8 °C.

Urines can be preserved with thymol to avoid bacterial action or by maintaining the pH below 4.

REFERENCE VALUES (1,2)

Sérum, plasma :	13 - 43	mg/dL
	2,1 - 7,1	mmol/L

Urine :	26 - 43	g/24 h
	0,43 - 0,71	mol/24 h
	1,7 - 2,9	g/dL*
	0,29 - 0,47	mol/L*

* for an urinary volume of 1.5 L per 24 hours

Note : It is recommended for each laboratory to establish and maintain its own reference values. The data given here are only an indication.

PROCEDURE

This reagent can be used on most analysers, semi-analysers and manual method.

The applications are available on request.

Wavelength : 340 nm

Temperature : 37 °C

Read against reagent blank.

- **One-reagent procedure**

	BLANC	CALIBRATION	TEST
Working Reagent	300 µL	300 µL	300 µL
Distilled water	3 µL	-	-
Standard	-	3 µL	-
Sample	-	-	3 µL

Mix and read the variation of absorbance (ΔA) between 25 seconds and 75 seconds.

- **Two-reagent procedure**

	BLANC	CALIBRATION</
--	-------	---------------

☛ In very rare cases, monoclonal gammopathies (multiple myeloma), in particular IgM type (Waldenström's macroglobulinemia) can cause unreliable results.⁽⁸⁾

Other compounds may interfere. ⁽⁹⁻¹⁰⁾

- On board stability / calibration frequency on Cobas Mira (no refrigerated)

On-board stability : 5 days for one-reagent procedure and 14 days for two-reagent procedure (capped vials and stored at 2-8 °C during the night)

Calibration frequency : 5 days for one-reagent procedure and 14 days for two-reagent procedure.

Recalibrate when reagent lots change, when quality control results fall outside the established range, and after a maintenance operation.

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- Newman, D. J., Price C. P., *Non protein Nitrogen Metabolite. Tietz Fundamentals of Clinical Chemistry*, 5th Ed., Burtis, C.A. & Ashwood, E.R. (W.B. Saunders eds. Philadelphia USA), (2001), 414.
- Tietz, N.W., *Clinical guide to laboratory tests*, 3rd Ed., (W.B. Saunders eds. Philadelphia USA), (1995), 622.
- First, M.R., *Renal function. Clinical Chemistry: Theory, Analysis, Correlation*, 4th Ed., Kaplan, L.A, Pesce, A.J., Kazmierczak, S.C., (Mosby Inc. eds St Louis USA), (2003), 477 and appendix.
- Bretaudière, J.P., et al., *Direct Enzymatic Determination of Urea in Plasma and Urine with a Centrifugal Analyzer. Clin. Chem.*, (1976), **22**, 1614.
- Fawcett, J.K., Scott, J.E., *A Rapid and Precise Method for the Determination of Urea. J. Clin. Path.*, (1960), **13**, 156.
- Vassault, A., et al., *Ann. Biol. Clin.*, (1986), **44**, 686.
- Vassault A., et al., *Ann. Biol. Clin.*, (1999), **57**, 685.
- Berth, M. & Delanghe, J., *Protein precipitation as a possible important pitfall in the clinical chemistry analysis of blood samples containing monoclonal immunoglobulins: 2 case reports and a review of literature, Acta Clin Belg.*, (2004), **59**, 263.
- Young, D.S., *Effects of preanalytical variables on clinical laboratory tests*, 2nd Ed., AACC Press, (1997).
- Young, D.S., *Effects of drugs on clinical laboratory tests*, 4th Ed., AACC Press, (1995).

SYMBOLS UTILISES SUR LES ETIQUETTES
 SYMBOLS USED ON LABELS
 SIMBOLOS USADOS EN LA ETIQUETA

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© : Modification par rapport à la version précédente/Modification from previous version