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AUML - 0420	6 x 50 mL
AUML - 0500	6 x 100 mL
AUML - 0700	4 x 250 mL
AUML - 0427	6 x 50 mL
AUML - 0507	6 x 100 mL
AUML - 0707	4 x 250 mL

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English - GB

In vitro diagnostic reagent, for professional use only

CLINICAL SIGNIFICANCE (1-3)


Uric acid is the major product of the catabolism of endogenous and exogenous (dietary) purine nucleosides (adenosine and guanosine). This transformation mainly occurs in the liver. Approximately 75% of uric acid is eliminated by kidneys; the remainder is secreted into the gastrointestinal tract, where it is degraded by bacterial enzymes. Uric acid is not very soluble in water; urate crystals can occur in urines when the concentration is abnormally high. It can also happen in plasma, crystals then deposit essentially in joints, which induce intense inflammatory responses (gout). Some causes for increasing uric acid rate in serum are: increasing of purines synthesis, metabolic disorders (Lesch-Nyhan syndrome for example), nutritional troubles, increasing of nucleic acid turn-over in case of proliferation of tumor cells, leukaemia, psoriasis, cytotoxic drugs, renal failures... Decreasing of uric acid rate in serum is more uncommon. It can occur in different cases: failure in renal elimination of uric acid (Fanconi syndrome), Hodgkin's disease for example.

METHOD (4)

Enzymatic - colorimetric. Trinder. End point.

PRINCIPE (4)

Enzymatic determination of uric acid according to the following reactions:

EHSPT = N-Ethyl-N-(2-Hydroxy-3-Sulfo-propyl) m-Toluidine
4-AAP = Amino-4-antipyrine

REAGENTS COMPOSITION

Reagent : R	
Phosphate buffer, pH 7.0	100 mmol/L
EHSPT	0,72 mmol/L
Ferrocyanide	0,03 mmol/L
Amino-4-antipyrine	0,37 mmol/L
Uricase	150 U/L
Peroxydase	≥ 12 000 U/L
Standard: Std	
Uric acid	6 mg/dL
	357 µmol/L

This standard is included in kits AUML-0427/0507/0707 and can be sold separately under the reference ACUR-0055.

MATERIAL REQUIRED BUT NOT PROVIDED

CONT-0060	ELITROL I	10 x 5 mL
CONT-0160	ELITROL II	10 x 5 mL

PRECAUTIONS

- Uric Acid Standard contains 0.15 % of sodium azide. It is harmful (Xn).

R22 : Harmful if swallowed.

S46 : If swallowed, seek medical advice immediately and show this container or label.

- Reagent contains less than 0.1 % sodium azide. Sodium azide can react with copper and lead plumbing to form explosive metal azides. If discharge in the canalisations, rinse with plenty of water.

- Use clean or single use laboratory equipment only to avoid contaminations.

- The standard should be immediately and tightly capped to prevent contamination and evaporation.

- For more information, Material Safety Data Sheet (MSDS) is available on request for professional user.

WASTE MANAGEMENT

Disposal of all waste material should be in accordance with local and legal requirements.

STABILITY OF REAGENTS

Stores at 2-8 °C and protect from light.

The reagent and the standard are stable until the expiry date stated on the label.

On board stability :

The stability is specific for each analyser (for COBAS MIRA refer to § PERFORMANCE DATA)

PREPARATION AND STABILITY OF WORKING REAGENT

The reagent and the standard are ready to use.

SAMPLES (1,3,5)

- Specimen

Sérum.
Heparinized plasma.

Urine diluted 1/10 in distilled water. If unpreserved urine is received, add 0.1 mL of 12.5M NaOH to 10 mL of well-mixed urine; mix well. Warming at 60 °C to dissolve precipitate may be needed.

- Storage

Sérum and heparinized plasma samples are stable 3 to 5 days if stored at 4 °C, 6 months at -20 °C. Urines are stable 3 days at room temperature. Do not refrigerate urine samples.

REFERENCE VALUES (3)

Man		Woman	
Sérum, plasma :	3,5 - 7,2	2,6 - 6,0	mg/dL
	208 - 428	155 - 357	µmol/L

Urine :	250 - 750	mg/24 h
	1,48 - 4,43	mmol/24 h
	16,7 - 50,0	mg/dL*
	0,99 - 2,97	mmol/L*

* for an urinary volume of 1.5 L per 24 hours

Note : It is recommended for each laboratory to establish and maintain its own reference values. The data given here are only an indication.

PROCEDURE

This reagent can be used on most analysers, semi-analysers and manual method.

The applications are available on request.

Wavelength : 550 nm

Temperature : 37 °C

Read against reagent blank.

	BLANK	CALIBRATION	TEST
Reagent R	200 µL	200 µL	200 µL
Distilled water	5 µL	-	-
Standard	-	5 µL	-
Sample	-	-	5 µL

Mix and measure the absorbances (A) after 325 second incubation.

CALCULATION

A Sample x n n = standard concentration
A Standard

Conversion factor : mg/dL x 59,48 = µmol/L
mg/dL x 0,059 = mmol/L
mg/dL x 10 = mg/L

Take the dilution factor into account for the calculation of uric acid concentration in urine.

CALIBRATION

Concentration value of Uric Acid Standard 6 mg/dL is traceable to the Standard Reference Material SRM 909b (of the National Institute of Standards and Technology).

The calibration frequency is specific for each analyser (for Cobas Mira refer to § PERFORMANCE DATA).

QUALITY CONTROL

To ensure adequate quality, control sera such as ELITROL I (normal control) and ELITROL II (abnormal control) are recommended.

PERFORMANCE DATA at 37 °C on COBAS MIRA

- Analytical range

The reagent is linear from 0.5 to 25 mg/dL (30 to 1487 µmol/L).

- Detection limit (6)

Determined according to SFBC protocol, the detection limit is equal to 0.2 mg/dL (12 µmol/L).

- Analytical sensitivity

The average variation of the analytical signal is 21x10⁻³ ΔA per mg/dL of uric acid (2.1x10⁻³ ΔA per mg/L, 350 ΔA per µmol/L) for a light path of 1 cm.

- Precision

Within-run reproducibility	n	Mean		CV (%)
		mg/dL	µmol/L	
Medium level	10	4,5	268	1,0
High level	10	10,6	631	0,7

Between-run reproducibility	n	Mean		CV (%)
		mg/dL	µmol/L	
Medium level	18	4,8	286	1,8
High level	18	8,9	529	2,7

- Correlation

A comparative study has been performed between Elitech method and another commercial reagent (enzymatic colorimetric -Uricase/PAP method) on 39 human serum samples. The sample concentrations were between 1.9 and 17.8 mg/dL (113 and 1059 µmol/L).

The parameters of linear regression are as follows :

Correlation coefficient : (r) = 0,9992

Linear regression : y = 0,9752 x - 0,24 mg/dL (14 µmol/L)

- Interferences (6)

According to SFBC recommendations, studies have been performed to determine the level of interference from different compounds :

Bilirubin : No significant interference up to 30 mg/dL (513,1 µmol/L).

Haemoglobin : Positive bias from 50 mg/dL (0,5 g/L).

Glucose : No significant interference up to 500 mg/dL (27,75 mmol/L).

Turbidity : Positive bias from 370 mg/dL (4,18 mmol/L) Triglyceride equivalent.

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